

ARGO-T (Transseptal)

CLOSED VIAL DISPENSING SYSTEM



- *DISPENSING IN CLOSED VIALS*
- *COMPLETE FILLING CYCLE IN 30 SECONDS*
- *PRECISION IN THE DISPENSED VOLUME*
- *GMP COMPLIANT*
- *DISPOSABLE STERILE KIT*

ARGO-T (TRANSSEPTAL) is an automatic dispensing system for closed vials able to fill a vial through rubber stopper. The equipment is designed to dispense radiotracers used in PET and SPECT diagnosis and metabolic radiotherapy.

ARGO-T is a system that complies with the cGMP guidelines for filling closed vials.

This device is part of the ARGO dispenser family with the unique feature to allow automatic perforation of the vial with a needle. Flexible holder support guarantee the possibility to easy adapt the dispenser with different vial format.

It is usually integrated in Class A shielded isolators equipped with manipulators, thus ensuring the utmost result in ergonomics and convenience in total compliance with the cGMP guidelines.

Main technical features

- Controlled volume dispensing by means of an integrated pumping syringe in the disposable kit
- Specific designed needle for piecing and filling together with integrated vent needle provided with 0.22 sterile filter
- Compact dimensions and low weight for easy moving and positioning in the work area
- Accurate volume dispensing within 5% for volumes from 1ml and up
- Label printer
- Barcode reader
- Calculation of the concentration and formulation of the initial volume of radiopharmaceuticals with the formulation unit
- Sterilising filter and bubble point test in-line.

Dati tecnici

| | ARGO TRANSSEPTAL |
|---------------------------------|---|
| Frame support material | AISI 316L - Mirror-Bright |
| External structure material | AISI 316L |
| External dimensions | 225 x206,5 x505 mm (w x d x h) |
| Weight | 15 kg |
| Volumetric dispensing accuracy | 0.5 ÷ 1 ± 0.05 ml 1 ÷ 10 ± 5% ml |
| Technical gas supply (He / N) | 5 bar |
| Type of containers to be filled | Closed vials |
| Dispensing method | Transseptal dispensing through rubber stopper |
| Minimum volume dispensed | 0.5 ml |

Kontakt / Contact:



COMECER S.p.A. - Via Maestri del Lavoro, 90
48014 - Castel Bolognese (RA) - Italy
t: +39 0546 656375 - f: +39 0546 656353
comecer@comecer.com - www.comecer.com

Part of COMECERGROUP



WINKGEN MEDICAL SYSTEMS GmbH & Co. KG
Eichendorffstraße 3
D-35584 Wetzlar
Telefon: +49(0)6441 / 381437
Telefax: +49(0)6441 / 381442
www.winkgen.de / email@winkgen.de



www.comecer.com/ARGO-T